

miglionico

DENTAL EQUIPMENT



LANGUAGE





ERGONOMICS
AND FUNCTIONALITY
AT YOUR FINGERTIPS

CONFIGURATION FLEXIBILITY

The NIL unit features a dentist's console on a double-jointed arm with 180° rotation.


The handle placed on the front of the console allows an easy grab-and-move movement.

The unit has a wide range of motion and is perfect to work single or four-handed.

The console can accommodate up to 5 instruments and it is equipped with an autoclavable silicone mat where to place a trayholder.

The capacitive touch control panel manages all the functions of the unit and the chair.

The instruments are placed in removable and autoclavable instrument holders with automatic activation sensors.



Configuration Flexibility



BASIC CONFIGURATION

CONFIGURATION
HYDRO CASE +
LIGHT POLE





**CONFIGURATION
WITH LIGHT POLE**



**CONFIGURATION
HYDRO CASE +
LIGHT POLE +
ASSISTANT ELEMENT**





1. COLD WATER CUP FILLER AND CUSPIDOR FLUSH ACTIVATION

4. INSTRUMENT WATER/SPRAY ACTIVATION

7. RECALL TRENDELENBURG POSITION

10. CHAIR MOVEMENTS

2. WARM WATER CUP FILLER AND CUSPIDOR FLUSH ACTIVATION

5. ASSISTANT CALL / DOOR OPENER

8. CUSPIDOR FLUSHING ACTIVATION

11. MICROMOTOR ROTATION REVERSAL

3. RINSING POSITION

6. MEMORIZED POSITIONS

9. RECALL MEMORIES

12. OPERATING LIGHT

Adaptive backrest

The backrest with motorized vertical excursion is available as an optional, and it can be managed by the control panel to ensure the correct positioning of the patient's head on the headrest, whether **pediatric** or **adult**.



Chair

The chair is designed to offer comfort and good working ergonomics. It is covered in skai medical leather for easy sanitization. The synchronized movement of the backrest and the seat ensures optimal and prevents the sliding effect. The anatomic backrest, with lumbar support, is available in 2 versions: standard for all working ergonomics and wide in memory foam for surgery or greater comfort of the patient, ideal for long-lasting treatments.




STANDARD BACKREST



XL MEMORY FOAM BACKREST

Colours

STANDARD COLOURS



SPECIAL COLOURS



Foot controls

The foot control is solid and well-balanced with a steel base to ensure stability while working. It is equipped with a retractable handle which makes it easy to move. It is available in the standard or push version, both with cable or wireless.





miglionico
DENTAL EQUIPMENT

Miglionico S.r.l.

Via Molise, Lotti 67/68 Z.I.

70021 - Acquaviva delle Fonti (BA) ITALY

P. Iva: 05306940726

Tel +39 080 759552

www.miglionico.eu

info@miglionico.net



Italy:

commerciale@miglionico.net

Eu & Extra Eu

export@miglionico.net

